



# STaR

## IHE Integration Statement

Product Version: **1.86**  
Revision: **3**  
Date : **September 2009**

ETIAM Head Office in Europe:  
ETIAM Technopole Atalante-Champeaux  
2, rue du Pierre-Joseph Colin  
35000 Rennes  
France

Phone: +33.(0)2.99.14.33.88  
Fax: +33.(0)2.99.14.33.80

E-mail: [info@etiam.com](mailto:info@etiam.com)  
[www.etiam.com](http://www.etiam.com)

ETIAM Head Office in America:  
ETIAM Corp.  
185 Alewife Brook Parkway, Suite 410  
Cambridge, MA 02138  
USA

Toll Free: +1-(877) 384-2662 (USA & Canada)  
Phone/Fax: +1-(617) 395-5809

<b>IHE Integration Statement</b>			
<b>Vendor</b>	<b>Product Name</b>	<b>Version</b>	<b>Date</b>
ETIAM	STaR	1.86	September 2009
This product implements all transactions required in the IHE Technical Framework to support the IHE Integration Profiles, Actors and Options listed below:			
<b>Integration Profiles Implemented</b>		<b>Actors Implemented</b>	<b>Options Implemented</b>
<b>Consistent Time</b>		Time Client	None
<b>Scheduled Workflow</b>		Patient Demographics Consumer	None
		Performed Procedure Step Manager	None
<b>Patient Information Reconciliation</b>		Image Manager/Image Archive	None
		Performed Procedure Step Manager	None
<b>Access to Radiology Information</b>		Image Manager/Image Archive	None
		Report Repository	None
<b>Key Image Note</b>		Image Manager/Image Archive	None
<b>Simple Image and Numeric Report</b>		Report Repository	None
<b>Consistent Presentation of Images</b>		Image Manager/Image Archive	None
<b>Reporting Workflow</b>		Image Manager/Image Archive	None
<b>Evidence Documents</b>		Image Manager/Image Archive	None
<b>Cross-enterprise Document Sharing for Imaging</b>		Imaging Document Source	None
<b>Retrieve Information Display</b>		Information Source	None
<b>Internet address for vendor's IHE information:</b> <a href="http://www.etiam.com/ihe">http://www.etiam.com/ihe</a>			
Links to Standards Conformance Statements for the Implementation			
DICOM	<a href="http://www.etiam.com/ihe">http://www.etiam.com/ihe</a>		
HL7			
Links to General Information on IHE			
In North America: <a href="http://www.ihe.net">http://www.ihe.net</a>	In Europe: <a href="http://www.ihe-europe.org">http://www.ihe-europe.org</a>	In Japan: <a href="http://www.jira-net.or.jp/ihe-j">http://www.jira-net.or.jp/ihe-j</a>	