

The logo for etiam, featuring the word "etiam" in white lowercase letters on a blue square background. The letter "i" has a small white graphic element above it.

etiam

Secure Medical Networking

ETIAM DICOM Izer

IHE Integration Statement



Product Version: 4.04
Revision: 8
Date : April 2013

ETIAM Head Office in Europe:
ETIAM Technopole S.A.S.U.
2, rue du Pierre-Joseph Colin
35000 Rennes
France

Phone: +33 (0)2 99 14 33 88
Fax: +33 (0)2 99 14 33 80

E-mail: info@etiam.com
www.etiam.com

ETIAM Head Office in America:
ETIAM Corp.
185 Alewife Brook Parkway, Suite 410
Cambridge, MA 02138
USA

Toll Free: +1-(877) 384-2662 (USA &
Canada)
Phone/Fax: +1-(617) 395-5809



IHE Integration Statement			
Vendor	Product Name	Version	Date
ETIAM	DICOM Izer	4.04	April 2013
This product implements all transactions required in the IHE Technical Framework to support the IHE Integration Profiles, Actors and Options listed below:			
Integration Profiles Implemented	Actors Implemented	Options Implemented	
Patient Demographics Query (PDQ)	Patient Demographics Consumer	Patient Demographics and Visit Query	
Portable Data for Imaging (PDI)	Media Creator	---	
Internet address for vendor's IHE information: http://www.etiam.com/ihe			
Links to Standards Conformance Statements for the Implementation:			
DICOM	http://www.etiam.com/ihe		
HL7			
Links to General Information on IHE			
In North America: http://www.ihe.net	In Europe: http://www.ihe-europe.org	In Japan: http://www.jira-net.or.jp/ihe-j	