



# ETIAM Viewer

## IHE Integration Statement

Product Version: **3.66**  
Revision: **7**  
Date: **June 2012**

ETIAM Head Office in Europe:  
ETIAM Technopole Atalante-Champeaux  
2, rue du Pierre-Joseph Colin  
35000 Rennes  
France

Phone: +33.(0)2.99.14.33.88  
Fax: +33.(0)2.99.14.33.80  
E-mail: [info@etiam.com](mailto:info@etiam.com)  
[www.etiam.com](http://www.etiam.com)

ETIAM Head Office in America:  
ETIAM Corp.  
185 Alewife Brook Parkway  
Suite 410  
Cambridge, MA 02138  
USA

Toll Free: +1-877-384-2662 (USA & Canada)  
Phone/Fax: +1-617-395-5809  
E-mail: [info@etiam.com](mailto:info@etiam.com)  
[www.etiam.com](http://www.etiam.com)

<b>IHE Integration Statement</b>			
<b>Vendor</b>	<b>Product Name</b>	<b>Version</b>	<b>Date</b>
ETIAM	ETIAM Viewer	3.66	June 2012
This product implements all transactions required in the IHE Technical Framework to support the IHE Integration Profiles, Actors and Options listed below:			
<b>Integration Profiles Implemented</b>	<b>Actors Implemented</b>	<b>Options Implemented</b>	
<b>Scheduled Workflow (SWF)</b>	Image Display	---	
<b>Portable Data for Imaging (PDI)</b>	Portable Media Importer	---	
	Image Display	---	
	Report Reader	---	
	Display	---	
Media Creator	---		
<b>Patient Demographics Query (PDQ)</b>	Patient Demographics Consumer	Patient Demographics and Visit Query	
<b>Import Reconciliation Workflow (IRWF)</b>	Importer	Unscheduled Import	
Internet address for vendor's IHE information: <a href="http://www.etiam.com/ihe">http://www.etiam.com/ihe</a>			
Links to Standards Conformance Statements for the Implementation:			
DICOM	<a href="http://www.etiam.com/ihe">http://www.etiam.com/ihe</a>		
HL7			
Links to General Information on IHE			
In North America: <a href="http://www.ihe.net">http://www.ihe.net</a>	In Europe: <a href="http://www.ihe-europe.org">http://www.ihe-europe.org</a>	In Japan: <a href="http://www.jira-net.or.jp/ihe-j">http://www.jira-net.or.jp/ihe-j</a>	