



Secure Medical Networking

ETIAM MARS

IHE Integration Statement



Product Version: 4.20
Revision: 9
Date: June 2013

ETIAM Head Office in Europe:
ETIAM S.A.S.U.
2, rue du Pierre-Joseph Colin
35000 Rennes
France

Phone: +33.(0)2.99.14.33.88
Fax: +33.(0)2.99.14.33.80

E-mail: info@etiam.com
www.etiam.com

ETIAM Head Office in America:
ETIAM Corp.
185 Alewife Brook Parkway, Suite 410
Cambridge, MA 02138
USA

Toll Free: +1-(877) 384-2662 (USA &
Canada)
Phone/Fax: +1-(617) 395-5809

IHE Integration Statement			
Vendor	Product Name	Version	Date
ETIAM	ETIAM MARS	4.20	June 2013
This product implements all transactions required in the IHE Technical Framework to support the IHE Integration Profiles, Actors and Options listed below:			
Integration Profiles Implemented	Actors Implemented	Options Implemented	
Portable Data for Imaging (PDI)	Media Creator	Web content	
Key Image Note (KIN)	Image Display*	None	
Internet address for vendor's IHE information:			
http://www.etiam.com/ihe			
Links to Standards Conformance Statements for the Implementation:			
DICOM	http://www.etiam.com/ihe		
HL7			
Links to General Information on IHE			
In North America:	In Europe:	In Japan:	
http://www.ihe.net	http://www.ihe-europe.org	http://www.jira-net.or.jp/ihe-j	

* The LiteBox viewer automatically recorded on the CD/DVD displays Key Image Notes.